Niton Europe GmbH · Joseph-Dollinger-Bogen 9 · 80807 München · Germany

## CE Declaration of Conformity

**Manufacturer Name:** Thermo NITON Analyzers LLC  
**Manufacturers Address:** 900 Middlesex Turnpike  
Billerica, Massachusetts, MA 01824, USA

**European Representative:** Dr. Björn Klaue, General Manager  
Niton Europe GmbH  
Joseph-Dollinger-Bogen 9  
80807 München, Germany

We declare that the product:

**Product Type:** Thermo Scientific Niton Analyzer  
**Model Number:** XL2  
**Product Options:** All

Has been tested by Parker Chomerics Test Services and is in conformance with the following standards and requirements:

### **EMC:**

EMC Directive:	2004/108/EC
EN61326-1: 2005	Electrical Equipment for Measurement, Control and Laboratory use, EMC Requirements
EN 61000-3-2 : 2006	Harmonic Current Emissions
EN 61000-3-3:1995 +A1:2001+A2: 2005	Voltage Fluctuations and Flicker
EN 61000-4-2: 1995+A1: 1998+A2: 2002	Electrostatic Discharge Immunity
EN 61000-4-3: 2002+A1: 2002	Radiated Electromagnetic Field Immunity
EN 61000-4-4 : 2004	Electrical Fast Transient Burst Immunity
EN 61000-4-5 : 2005	Surge Immunity Requirements
EN 61000-4-6 :	Radio Frequency Common Mode Immunity
EN 61000-4-11 : 2004	Voltage Dips, Short Interruptions and Voltage Variations
CISPR 11:	Limits and Methods of Measurements of Radio Disturbance Characteristics of Industr., Scient. and Medical (ISM) Radio Frequency Equipment Group 1, Class A

**Test Laboratory:** Parker Chomerics Test Services, 77 Dragon Court, Woburn, Massachusetts 01801. Test report # EMI5408.09.

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
# Niton XL2 Analyzer Quick Start Guide

**IMPORTANT SAFETY INFORMATION:**  
We recommend that you read the entire user's manual before using your analyzer.

**1**


**Connect the Charger To the Wall Outlet**

**Orient the Battery Pack Properly**

MAKE SURE THAT THE KEY IS IN FRONT AND THE THIN SIDE IS DOWN

GREEN LIGHT INDICATES FULLY CHARGED


**Place Battery in Charger**

**2**

PUSH BATTERY PACK INTO HANDLE, THIN SIDE UP, THEN SECURE LATCH

**Insert Battery**


**Analyzer Controls**


PRESS AND HOLD ON/OFF/ESCAPE BUTTON UNTIL TOUCH SCREEN POWERS ON

**Power Up**


**3**

**Startup Screens**

**Perform System Check**


**Main Menu**


# Niton XL2 Analyzer Quick Start Guide


**IMPORTANT SAFETY INFORMATION:**  
We recommend that you read the entire user's manual before using your analyzer

**Select Sample Type For Analysis**


**From Ready To Test Screen, Analyze Alloy Check Sample Provided**


**Stop Reading After 10 Secs**


REFER TO THE USER'S GUIDE FOR SPECIFIC STANDARD OPERATING PROCEDURE


**Install PC Software Resource Guide/ User's Guide**


**Start PC Software**

**Install USB Driver Provided On Disk**

**Connect to PC**

**Download Data to PC**

**Data Report**  
PART NUMBER 410-007

Niton Europe GmbH · Joseph-Dollinger-Bogen 9 · 80807 München · Germany

## CE Declaration of Conformity

**Manufacturer Name:** Thermo NITON Analyzers LLC  
**Manufacturers Address:** 900 Middlesex Turnpike  
Billerica Massachusetts, 01824  
USA

**European Representative:** Björn Klaue, General Manager  
Niton Europe GmbH  
Joseph-Dollinger-Bogen 9  
80807 München

**Declares that the product:**

**Product Type:** ThermoFisher Scientific Niton Analyzer  
**Model Number:** XL3t  
**Product Options:** All

**Has been tested and is in conformance with the following standards and requirements:**

**EMC:**  
EN61326 Electrical Equipment for Measurement, Control and Laboratory use, EMC Requirements  
EN 61000-4-2 Electrostatic Discharge  
EN 61000-4-3 Radiated Electromagnetic Fields  
EN 61000-4-4 Electrical Fast Transient/Burst  
EN 61000-4-5 Surge Immunity Requirements  
EN 61000-4-6 Conducted Disturbances Induced by Radio Frequency Fields  
EN 61000-4-8 Power Frequency Magnetic Field Immunity Test  
EN 61000-4-11 Voltage Dips, Short Interruptions and Voltage Variations  
EN 55011 Limits and Methods of Measurements of Radio Disturbance Characteristics of Industrial, Scientific and Medical (ISM) Radio Frequency Equipment Group 1, Class A

**Safety:**  
IEC 61010-1 2001 (Second Edition)  
Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory use

München, November 7<sup>th</sup>, 2007


Dr. Björn Klaue  
General Manager  
Director Technical Support and Applications  
Radiation Safety Officer

**Niton Europe GmbH**  
Joseph-Dollinger-Bogen 9 · D-80807 München  
Tel. +49-89-36 61 38 0 · Fax +49-89-36 61 33 0  
E-mail: europe@niton.com

Niton Europe GmbH  
Joseph-Dollinger-Bogen 9  
80807 München  
Germany  
Tel. +49-(0)89-36 61 38-0  
Fax +49-(0)89-36 61 38-30  
Email niton\_eur@thermofisher.com

Bankverbindung  
Commerzbank AG  
Kto-Nr 8218000  
BLZ 763 400 61  
IBAN DE 58 7634 0051 0921 8000 00  
SWIFT COBADEFF763  
Geschäftsführer:  
Henry J. Grodzins, Dr. Olaf Haupt  
Dr. Björn Klaue, James R.E. Coley,  
Piet van der Zande  
Amtsgericht München HRB 129790  
Ust-ID-Nr.: DE 205 372 303

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


Version 6.0

**Thermo**  
SCIENTIFIC

# NITON XL3t Quick Start and Radiation Safety Guide

## Unpacking and Assembling Your NITON XRF Analyzer

- Inspect the shipping carton for signs of damage such as crushed or water damaged packaging. Immediately notify the shipping company and Thermo Fisher Scientific, in the United States, toll free, at (800) 875-1578, or outside the United States, at +1-978-670-7460, if there is any visible damage to the shipping container or any of its contents.
- Open the packing carton. If your analyzer is not packed in its carrying case, please call Thermo Fisher Scientific immediately, in the United States, toll free, at (800) 875-1578, or outside the United States, at +1-978-670-7460.
- Verify the contents of the shipping container against the enclosed packing list. If there are any discrepancies between the actual contents of the shipping container and the enclosed packing list, please notify Thermo Fisher Scientific immediately, in the United States, toll free, at (800) 875-1578, or outside the United States, at +1-978-670-7460.
- Open the carrying case and visually inspect the analyzer for damage before removing it from the case. Call the shipper and Thermo Fisher Scientific if you find any damage to the case or its contents.
- Save the shipping carton and all packing materials. Store them in a safe, dry area for reuse the next time that you ship the analyzer.

## The NITON XRF Analyzer Overview

The NITON XL3 Analyzer is a single unit, hand held, high performance portable x-ray fluorescence (XRF) elemental analyzer.



Figure 0-1. Analyzer Overview

## The Control Panel

The control panel is located on the analyzer's top housing, directly below the LCD touch screen (see Figure 0-1). The control panel consists of a 4 way touch pad, a center button, and two control buttons, one on each side. Using either the control panel or the touch screen you may navigate through all of the analyzer's screens and menus. You can control the movement of the screen cursor by pressing the four way control pad in one of four directions to highlight each of the menu options. The Select button in the center of the four way touch pad is used to select highlighted menu options. The on/off/escape button both controls the power to the analyzer and serves as an "escape" button. When the on/off/escape button is pushed and immediately released, it functions as an "escape", and brings you back to the Main Menu from the current screen in the menu system.

**Figure 0-2. The Control Panel**


Push and hold the on/off/escape button for at least 3 seconds to turn the analyzer on. Push the on/off/escape button and hold it down for about 10 seconds to shut off power to the analyzer from any screen in the menu system.

You also have the option of operating the analyzer, including navigating the menu system, by using the built in touch screen. To select a menu option, tap on the icon once. The touch screen icons have the same functionality as the four way touch pad, the on/off/escape button, and the select or enter button. This Quick Start will refer to the process of choosing a course of action by selecting an icon from a menu, either using the touch screen or using the control panel buttons, as “selecting.”

Selecting the **Return** icon works everywhere throughout the User Interface to bring you back to the previous menu from the current menu in the menu system. Use the on/off/escape button to return to the **Main Menu**.

## The LCD Touch Screen

The LCD Touch Screen on your XL3 Analyzer is designed to swing up and down to different angles for ease in viewing and interacting with your analyzer. The LCD Touch Screen is connected to your analyzer along the base of the screen, right above the Control panel. The screen is not designed to separate from the analyzer, but can be adjusted to any arbitrary angle between zero degrees - that is, flush with the analyzer - and 85 degrees, which is almost perpendicular. The LCD Touch Screen will stay at any given angle between these extremes until moved to a different angle. When in closed position, the screen is secured by a catch at the top center of the screen housing.

**Figure 0-3. XL3 Analyzer Showing LCD Screen Tilted.**

- To raise the LCD Touch Screen, disengage the catch at the top-center of the LCD Touch Screen housing and gently pull the screen towards you until it is at the best angle for your use.
- To close the LCD Touch Screen, gently push away from you along the top edge of the screen housing. The screen will swing down until the catch solidly engages with an audible click.

**Note** The LCD Touch Screen cannot be removed from your XL3 analyzer. Removing or attempting to remove the LCD Touch Screen will damage your analyzer and void your warranty.

**Note** Always close your LCD Touch Screen before storing or transporting your XL3 analyzer.

## The Battery Life Indicator

The Battery Life Indicator is visible on all screens in the menu system. The indicator is visible in the top right portion of the screen, and graphically shows you how much battery life is left, enabling you to change batteries as needed to avoid unexpected shutdowns.


Figure 0-4. Battery Life Indicator

The more green visible in the indicator, the higher the charge. The more red visible in the indicator, the lower the charge. It's best to charge one battery while using the other, to avoid work slowdowns or stoppages due to battery charge conditions.

## Radiation and General Safety

This chapter covers topics related to radiation safety and general safety when using a Thermo Scientific NITON XL3t analyzer. At a minimum all operators of the XL3t should be familiar with the instructions provided in this chapter in order to handle the XL3t in a safe manner. In addition to reading the information presented on the following pages, Thermo Fisher Scientific recommends that instrument users participate in a radiation safety and operational training class.

**WARNING!** Always treat radiation with respect. Do not hold your analyzer near the Kapton window during testing. Never point your analyzer at yourself or anyone else when the shutter is open. ♦

## Radiation Protection Basics

The NITON Model XL3t analyzer contains an x-ray tube which emits radiation only when the user turns the x-ray tube on. When the x-ray tube is on and the shutter is open, as during a measurement, the analyzer emits a directed radiation beam (See Figures 0-6 and 0-7). Reasonable effort should be made to maintain exposures to radiation as far below dose limits as is practical. This is known as the ALARA (As Low as Reasonably Achievable) principle. For any given source of radiation, three factors will help minimize your radiation exposure: Time, Distance, and Shielding.

**Time** The longer you are exposed to a source of radiation the longer the radiation is able to interact in your body and the greater the dose you receive. Dose increases in direct proportion to length of exposure.

**Distance** The closer you are to a source of radiation, the more radiation strikes you. Based on geometry alone, dose increases and decreases with an inverse-squared relation to your distance from the source of radiation (additional dose rate reduction comes from air attenuation). For example, the radiation dose one foot from a source is nine times greater than the dose three feet from the source. Remember to keep your hands and all body parts away from the front end of the analyzer when the shutter is open to minimize your exposure.

**Shielding** Shielding is any material that is placed between you and the radiation source. The more material between you and the source, or the denser the material, the less you will be exposed to that radiation. Supplied or optional test stands are an additional source of shielding for analysis. A backscatter shield accessory is also available and may be appropriate in some applications.

## Exposure to Radiation

Human dose to radiation is typically measured in rem, or in one-thousandths of a rem, called millirem (mrem), 1 rem = 1000 mrem. Another unit of dose is the Sievert (Sv), 1 Sv = 100 rem. The allowable limit for occupational exposure in the U.S (and many countries internationally) is 5,000 mrem/year (50 mSv/year) for deep (penetrating) dose and 50,000 mrem/year (500 mSv/year) for shallow (i.e., skin) dose or dose to extremities. Deep, shallow, and extremity exposure from a properly used NITON XL3t analyzer should be less than 200 mrem per year, (2.0 mSv per year) even if the analyzer is used as much as 2,000 hours per year, with the shutter open continuously. The only anticipated exceptions to the 200 mrem maximum annual dose are: 1) routine and frequent analysis of plastic samples without use of a test stand, backscatter shield, or similar additional protective measures, or 2) improper use where a part of the body is in the primary beam path. **NEVER OPERATE THE DEVICE WITH A PART OF YOUR BODY IN THE PRIMARY BEAM PATH OR WITH THE PRIMARY BEAM PATH DIRECTED AT ANYONE ELSE.** Also, consider the use of protective accessories such as a shielded test stand or backscatter shield (or equivalent) when performing routine and/or frequent analysis of any of the following:

- plastic (or similarly low density) samples,
- thin samples (such as foils, circuit boards, and wires), or
- samples that are smaller than the analysis window.

Shown in Table 0-1 below are the typical background radiation doses received by the average member of the public. The radiation dose limits for radiation workers in are also shown in Table 0-2.

## How to Use the NITON XL3t Analyzer Safely

The NITON XL3t analyzer is designed to be safe to operate provided that it is used in accordance with manufacturers' instructions. Under conditions of normal use, monitored operators seldom receive a measurable dose and have not been known to receive in excess of 10% of the annual occupational dose limits (a criteria that would require monitoring under regulation in the U.S.). In addition to proper use of the XL3t, it is recommended that you follow these precautions to ensure your safety and the safety of those around you.

### Know where the beam is

The primary beam is a directed beam out of the front of the analyzer that can have high dose rates. The secondary beam, or scattered beam, has much lower dose rates.



Figure 0-6. . Primary Beam


**Figure 0-7. Secondary (Scattered) Beam**

### **The Shutter-Open Indicator Lights**

When the lights are flashing, the primary beam is on, the shutter is open, and radiation is being emitted from the front of the analyzer. (This does not include the brief flash of the lights when first turning the analyzer on.)

**Figure 0-8. X-ray Beam Indicator Lights**

**Handle and Use with Respect**

Avoid holding the front of the analyzer when the x-ray tube is energized and the shutter is open. Never point the instrument at yourself or anyone else when the shutter is open and the x-ray tube is energized. Never look into the path of the primary beam.

**Follow a Radiation Protection Program**

Your organization should establish, document, and follow a Radiation Protection Program. An example of such a program can be found on the resource CD (provided with the instrument).

**Take Proper Care of your NITON XL3t Analyzer**

Keeping your analyzer maintained in good condition will help minimize the risk of accidental exposure. Mechanical malfunction of the shutter can be avoided by maintaining the Kapton window, as described in the User Guide. This prevents foreign objects from entering your XL3t.

**Avoid Over-Exposures**

Direct contact with the window could result in overexposures in the times indicated in Table 0-4 below.

**Table 0-1. Potential Exposure Limit Times**

Location of Dose	Limit	Time to Reach Limit
Deep Dose / Whole Body	5 rem (50 mSv)	2.1 minutes
Shallow Dose / Extremities	50 rem (500 mSv)	0.95 minutes
Member of Public (i.e. untrained operator)	0.1 to 5 rem (1 to 50 mSv)	2.5 to 9.5 seconds

Extremity is defined by the NRC as the hand, elbow, arm below the elbow, foot, knee, or leg below the knee. Whole Body is defined by the NRC as the head, trunk (including male gonads), arms above the elbow, or legs above the knee.

## Safe Handling of Samples

As mentioned many times in this chapter, never place any part of your body in the path of the x-ray beam. There is always a safe way to handle samples whether they are small, irregularly shaped, or of low density. Never look into the path of the primary beam.

### Small Samples

A small sample would be any sample that is smaller than the Kapton window. Small samples present a unique risk because they don't block the entire beam path. The difficulty with placing small samples down on a work surface to analyze them is that you may get readings from the work surface that interfere with analytical results. A test stand is an effective way of analyzing small samples accurately and safely. Never hold samples during analysis or look into the path of the primary beam.

### Irregularly Shaped Samples


Irregularly shaped samples may not allow the proximity button to be depressed, or they may not entirely cover the primary beam and cause additional scattering. A back scatter shield is a safe way of reducing your radiation exposure while effectively analyzing an irregularly shaped sample.

### Low Density Materials (such as plastics)


X-rays are attenuated more through denser materials and less through low density materials such as plastic. This causes higher dose rates in the scattered radiation. If you are frequently handling low density samples, you should consider the use of test stands, backscatter shields, or the equivalent.

# Niton XL2 Analyzer Quick Start Guide


**IMPORTANT SAFETY INFORMATION:**  
We recommend that you read the entire user's manual before using your analyzer

**Select Sample Type For Analysis**


**From Ready To Test Screen, Analyze Alloy Check Sample Provided**


**Stop Reading After 10 Secs**

REFER TO THE USER'S GUIDE FOR SPECIFIC STANDARD OPERATING PROCEDURE


## NDT DATA DOWNLOAD


**Install PC Software Resource Guide/ User's Guide**

**Start PC Software**


**Install USB Driver Provided On Disk**

**Connect to PC**

## QUERY READINGS


**Download Data to PC**

**Data Report**  
PART NUMBER 410-007


# Niton XL2 Analyzer Quick Start Guide

**IMPORTANT SAFETY INFORMATION:**  
We recommend that you read the entire user's manual before using your analyzer

**1**


**Connect the Charger To the Wall Outlet**

**Orient the Battery Pack Properly**

MAKE SURE THAT THE KEY IS IN FRONT AND THE THIN SIDE IS DOWN

**Place Battery in Charger**


GREEN LIGHT INDICATES FULLY CHARGED

**2**


PUSH BATTERY PACK INTO HANDLE, THIN SIDE UP, THEN SECURE LATCH

**Insert Battery**


**Analyzer Controls**


**Power Up**

PRESS AND HOLD ON/OFF/ESCAPE BUTTON UNTIL TOUCH SCREEN POWERS ON

**3**

**Startup Screens**

**Perform System Check**

**Main Menu**  
PART NUMBER 410-007